

FLUXANA

Your supplier for
XRF Application Solutions

HD ELEKTRONIK

Fusion machines

Liquid Analysis Technology



FLUXANA GmbH & Co. KG
Borschelstr. 3, 47551 Bedburg-Hau, Germany
Tel.: +49 (0) 2821 997 32-0
Fax: +49 (0) 2821 997 32 29
E-mail: info@fluxana.de
Web: www.fluxana.com

Amtsgericht Kleve: HR-A 2935, HR-B 8211
Ust-IdNr.: DE 814692564, Steuer-Nr. 116/5755/0442
Finanzamt Kleve

In cooperation with:

HD Elektronik und Elektrotechnik GmbH
Tichelstr. 10, 47533 Kleve, Germany
Tel.: +49 (0) 2821 148 10
Fax: +49 (0) 2821 148 09
E-mail: hde@hdelektronik.de
Web: www.hdelektronik.de

Amtsgericht Kleve: HR-B 2162
Ust-IdNr.: DE 812941185, Steuer-Nr. 116/57071869
Finanzamt Kleve

Official agent



NEW! Cups with film
Already prepared

Please ask for a quote.




Thin Films

The preferences of our thin films:

- 25% more pre-cut films per package than other brands
- dust proof anti static packaging
- lot number tracking system
- rolls, pre-cut circles & sheets



Item No.	Description	Thickness μm	Packaging
TF-025	Microporous		Roll 6,4cm x 5,1m
TF-112	Mylar®	12	Roll 7,6cm x 91,4m
TF-115	Ultrapolyester	1,5	Roll 7,6cm x 91,4m
TF-125	Mylar®	2,5	Roll 7,6cm x 91,4m
TF-125-345	Mylar®	2,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-12560	PolyRep® (replacement for polycarbonate)	Mylar® 2,5 + PP 6,0	2 Rolls a 7,6cm x 91,4m
TF-135	Mylar®	3,5	Roll 7,6cm x 91,4m
TF-135-255	Mylar®	3,5	Circle 6,4cm (500 pcs.)
TF-135-355	Mylar®	3,5	Circle 8,9cm (500 pcs.)
TF-135-345	Mylar®	3,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-160	Mylar®	6,0	Roll 7,6cm x 91,4m
TF-160-255	Mylar®	6,0	Circle 6,4cm (500 pcs.)
TF-160-355	Mylar®	6,0	Circle 8,9cm (500 pcs.)
TF-160-345	Mylar®	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-212	Polypropylene	12,0	Roll 7,6cm x 91,4m
TF-212-345	Polypropylene	12,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240	Polypropylene	4,0	Roll 7,6cm x 91,4m
TF-240-255	Polypropylene	4,0	Circle 6,4cm (500 pcs.)
TF-240-30	Polypropylene	4,0	Circle 7,6cm (100 pcs.)
TF-240-345	Polypropylene	4,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240-355	Polypropylene	4,0	Circle 8,9cm (500 pcs.)
TF-250	Polypropylene	5,0	Roll 7,6cm x 91,4m
TF-260	Polypropylene	6,0	Roll 7,6cm x 91,4m
TF-260-255	Polypropylene	6,0	Circle 6,4cm (500 pcs.)
TF-260-355	Polypropylene	6,0	Circle 8,9cm (500 pcs.)
TF-260-345	Polypropylene	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-350	Polycarbonate	5,0	Roll 7,6cm x 91,4m
TF-475	Kapton®	7,5	Roll 7,6cm x 91,4m
TF-475-345	Kapton®	7,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-500	XRF Tape	-	Roll 25mm x 66m
TF-LS-240-255	Polypropylene Low Sulfur	4,0	Circle 6,4cm (500 pcs.)
WZ-0001	Preparation tool for our sample cups	-	Measurement: 56 x 20mm



Professional Cups

The sample cups of FLUXANA are designed for use in XRF. All cups have been optimized and approved by our customers world wide.



We offer different XRF-sample cups to you.

Item No.	Description	OD mm	ID mm	H mm	Vol ml
SC-3332	Cup, ring, cap for: Spectro, XEPOS, XLAB 2000®	31	24	23	9
SC-3335	Cup, ring, cap for Philips/Panalytical®	35	28	35	17
SC-3340	Inner ring, outer ring, cap for: Bruker®, SRS, S4, S8, ARL Optimix®	41	34	25	17
SC-3340T	Inner ring, outer ring, cap for: Bruker® SRS, S4, S8	41	32	36	22
SC-3345	Inner ring, outer ring, cap universal for 45mm	45	37	37	33
SC-3351	Cup, cap universal for 51,5mm holder	51	44	41	48
SC-4031	Cup, ring, single open ended, universal open for 32mm holder	31	24	24	9
SC-4047	Cup, ring, single open ended, for Horiba®	46	35	18	12
SC-4131	Cup, ring, collar single open ended, universal for 32mm holder	32	25	24	9
SC-4140	Cup, ring, collar single open ended, universal for 40mm holder	40	32	23	16
SC-4231	Cup, 2 rings, double open ended, universal for 32mm holder	31	24	24	10
SC-4240	Cup, 2 rings, double open ended, universal for 40mm holder	39	32	23	18
SC-4331	Cup, ring, cap universal for 32mm holder	31	24	26	9
SC-4335	Cup, ring, cap for Philips/Panalytical®	34	27	36	17
SC-4340	Cup, ring, cap universal for 40mm holder	39	31	25	15
SC-4345	Cup, cap for Philips/Panalytical®	45	38	41	35
SC-6232	Cup, 2 rings, double open ended for Spectro XEPOS, XLAB 2000®	32	26	21	11
SC-6332	Cup, ring, cap for Spectro XEPOS, XLAB 2000®	32	26	23	11
SC-7332	Inner cup and cap, designed for Oxford®	34	28	38	16
SC-7332inner	Only inner cup, designed for Oxford LABX®	32	28	36	22
SC-7332cap	Only cap, designed for Oxford LABX®	-	-	-	-
SC-8047	Cup, ring, single open ended, for Horiba® 47mm holder	42	36	19	18
SC-MC-1520	XRF micro sample cups, double open ended with vented caps	31	15	27	3,6

One set contains 100 pieces. Caps and rings are packed separately! No cutting necessary!

Certified Reference Materials

Catalog	Description
	Petro standards liquid (ICP / XRF) This catalog contains: <ul style="list-style-type: none"> - sulfur in oil - sulfur in gasoline - lubricants - wear metals - metal concentrates - XRF oil package - XRF monitors
	Waterbased standards for AA and ICP This catalog gives an overview of liquid reference materials used for ICP / ICP-MS, IC and AA: <ul style="list-style-type: none"> - ICP / ICP-MS standards - EPA standards - Ion chromatography standards - Atomic absorption standards